

Paper Presentation

Rapid, artifact-reduced, image reconstruction for
super-resolution structured illumination microscopy

B10901151 林祐群

Instructor T. Tony Yang

Outline

- A. Overview
- B. Background - Structured Illumination Microscopy
 - From Wiener SIM to HiFi SIM (RL Algorithm)
 - JSFR-SIM
- C. JSFR-AR-SIM
- D. Results & Analysis
- E. Merits & Future Possibility
- F. References


Overview

Rapid, artifact-reduced, image reconstruction for super-resolution structured illumination microscopy


Wang Z, Zhao T, Cai Y, et al.

Published April 13, 2023

- Combining high-speed reconstruction framework with high fidelity optimization approach designed to suppress the **sidelobe** artifact.




Reference: <https://www.mdpi.com/2076-3417/6/11/359>



- Achieve a **real-time** (millisecond-level, not include parameters estimation) reconstruction, display, and storage of high-fidelity, super-resolved images of the samples is achieved

Demo of the proposed method

The real-time reconstruction and display results of Argo-SIM slide.



Argo-SIM Slides



Reference: <https://www.sciencedirect.com/science/article/pii/S266667582300053X>

Reference: <https://www.klv.co.jp/product/slide-microscope/argo-hm-v2.html>

Structured Illumination Microscopy

Merits:

- A kind of **Far-field** super-resolution fluorescence microscopy
 - => Larger Field of View, simpler operation
 - => good for to visualize the fine details of intracellular structures.
- Enables high-speed imaging while **minimizing photodamage**
 1. Lower Illumination Intensity
 2. Wide-Field Illumination

Current Drawbacks:

- Reconstruction Artifacts
- SIM image acquisition and reconstruction process are separated (online v.s. offline) and time-consuming




Wiener SIM

Parameters Estimation

Input Images * 9
(3 Frequency bands)

Spectrum Separation,
Translation, & Summation


Filtered and
Inverse Fourier Transform



Suppress noise!
Increase resolution!

Reference: <https://www.nature.com/articles/s41377-023-01204-4>

Wiener Filter





A kind of MMSE filter!

Reference: <https://www.youtube.com/watch?v=Nm0UXL1OpRE>

Sidelobe Artifacts


Smaller modulation factors will cause the high-order shifted components in the recombined spectrum to be over-amplified. The equivalent OTF of the Wiener-SIM protocol exhibits **residual peaks and downward kinks** at the middle frequency, which is in essence different from the ideal OTF with doubled resolution.



Reference: <https://www.nature.com/articles/s41377-021-00513-w#Sec15>

HiFi SIM

The difference between the HiFi- and Wiener-SIM is mainly in the **preprocessing** procedure (iterations set to 5) and final **two-step optimization $\tilde{W}_1(\mathbf{k})$ and $\tilde{W}_2(\mathbf{k})$** . It can achieve **~2-fold of the diffraction limit**



Reference: <https://www.nature.com/articles/s41377-023-01204-4>


Preprocessing - Richardson–Lucy (RL) deconvolution

An **iterative** procedure for recovering an underlying image that has been **blurred by a known PSF**

$$\circ \hat{I}^{(t+1)} = \hat{I}^{(t)} \cdot \left(\frac{B}{\hat{I}^{(t)} \otimes K} \otimes K^* \right)$$

- Increase** the intensity of the **high-frequency** component
- Reduce the out-of-focus signal-related artifacts
 - Enhancing detection of reconstruction parameters.

Take Iteration...



| | |
|-------------|--------------------------------|
| \hat{I} : | <i>Latent Image Estimation</i> |
| B : | <i>Blurred Image</i> |
| K : | <i>Point Spread Function</i> |
| K^* : | <i>Flip PSF</i> |
| \otimes : | <i>2D Convolution</i> |

Reference: <https://www.youtube.com/watch?v=v1kE3mbexko>


Reference: https://scikit-image.org/docs/0.23.x/auto_examples/filters/plot_deconvolution.html

Two-Step Optimization: $\tilde{W}_1(\mathbf{k})$, $\tilde{W}_2(\mathbf{k})$

The difference between the HiFi- and Wiener-SIM is mainly in the **preprocessing** procedure (iterations set to 5) and final **two-step optimization** $\tilde{W}_1(\mathbf{k})$ and $\tilde{W}_2(\mathbf{k})$

$\tilde{W}_1(\mathbf{k})$:
>> to correct $\tilde{H}_{Theoretical}(\mathbf{k})$ to be close to
 $\tilde{H}_{Ideal}(\mathbf{k})$

$\tilde{W}_2(\mathbf{k})$:
>> to recover the high-frequency signals suppressed by the optimized OTF



However, both preprocessing and optimization are of huge computational complexity!!!

Reference: <https://www.nature.com/articles/s41377-021-00513-w#Sec15>

JSFR SIM

It is a **spatial domain-based** SIM reconstruction instead of a Fourier domain-based method. The reconstruction time is **10.2 ms** for 512*512 raw images









Parameters Estimation

Input Images * 9
(3 Frequency bands)

Filtering

Calculating coefficient matrix &
multiplying filtered raw images

Summing (1-9) &
Wiener Filtering



**Super-resolution
image**

Reference: <https://www.nature.com/articles/s41377-023-01204-4>

JSFR AR SIM

- HiFi- and JSFR-SIM both employed the OTF compensation and attenuation to enhance the spatial resolution and suppress the background fluorescence.
 - => This is equivalent to applying a **bandpass filter** $a(\mathbf{k})\tilde{H}^*(\mathbf{k})$ on all the preprocessed raw images
 - => The extra spectrum information gained by RL deconvolution will be **suppressed** by the bandpass filter
 - => RL deconvolution **doesn't provide extra** spectrum information of the sample than Wiener-type filter

Replace the RL deconvolution with Wiener-type filter to accelerate!


$\tilde{H}(\mathbf{k})$: OTF of the microscope

- Combine the two optimization steps in HiFi-Sim into a single one, since the two optimization functions $\tilde{W}_1(\mathbf{k})$ and $\tilde{W}_2(\mathbf{k})$ in HiFi-SIM do not vary with the sample distribution.


$$\tilde{W}(\mathbf{k}) = \tilde{W}_1(\mathbf{k}) \cdot \tilde{W}_2(\mathbf{k})$$

- The parameters of the illumination field (wave vectors, initial phases, and modulation depth) are **pre-computed** for later use, along with coefficient functions $c_{d,i}(\mathbf{r})$, composite filter $\tilde{G}(\mathbf{k})$, and the combined optimization function $\tilde{W}(\mathbf{k})$.

JSFR AR SIM




Quality Results - Argo-SIM




Quality Results - Microtubules


D







F



E



G



Speed Results - Computation Time

| Input image size | Output image size | Acquisition time (ms) ^a | Execution time of JSFR-AR-SIM (ms) | | Execution time of HiFi-SIM (ms) | |
|------------------|-------------------|------------------------------------|------------------------------------|--------------------|---------------------------------|----------------|
| | | | CPU ^b | GPU ^c | CPU | GPU |
| 2,048 × 2,048 | 4,096 × 4,096 | 95.0 | 5,940.1 ± 776.1 (9.9) ^d | 88.7 ± 12.3 (13.7) | 58,806.2 ± 1,290.8 | 1,205.6 ± 13.5 |
| 1,024 × 1,024 | 2,048 × 2,048 | 45.0 | 1,340.9 ± 70.8 (10.2) ^d | 21.5 ± 6.1 (16.4) | 13,720.8 ± 616.0 | 352.2 ± 6.3 |
| 512 × 512 | 1,024 × 1,024 | 22.5 | 274.9 ± 7.3 (12.4) | 5.7 ± 4.5 (17.3) | 3,412.6 ± 186.6 | 98.5 ± 3.7 |
| 256 × 256 | 512 × 512 | 11.3 | 68.1 ± 1.7 (12.7) | 3.1 ± 0.1 (10.0) | 866.2 ± 14.3 | 31.1 ± 0.6 |

JSFR-AR-SIM reconstructs SR images 2- to 4-fold more rapidly than it takes to acquire them

1. Platform: MATLAB
2. Personal computer (i7-9700K@3.6GHz, NVIDIA GeForce RTX 3080ti 12GB, Samsung 860 EVO 500GB SSD)

Results & Analysis - Workflow

HiFi-SIM

2.13s

Calculate attenuated OTF & Weight Functions

10.37s (~78%)

R-L deconvolution for each raw image

1.51s

Separate components in each pattern orientation

0.54s

Multiply each component with attenuated OTF

0.86s

Shifts and combines all the components

0.01s

Multiply the spectrum with the weight function

Inverse Fourier Transform

HiFi-SIM: 292*FT + 241*PM + 30*PA + 6*PP

JSFR-AR-SIM: 20*FT + 19*PM + 8*PA

Estimation the illumination parameters

JSFR-AR-SIM

Calculate Coefficient Functions, Composite Filter, Weight Functions

Preprocessing: background subtraction, upsampling, etc.

Filter all raw images with composite filter

Multiply with coefficient functions and sum all the product

Transform to frequency domain and multiply with weight functions

Inverse Fourier Transform

FT: (Inverse Fourier Transform) PM: Pointwise Multiplication

PA: Pointwise Addition

PP: Pointwise Power

Reference: <https://www.cell.com/cms/10.1101/j.xinn.2023.100425/attachment/2ea5b0e0-1a28-4360-8b50-26c2ee0fb6e6/mmc1.pdf>

Merits & Future Possibility

- High Speed + High-fidelity in a single framework.
=> achieve **real-time** reconstruction with **no sidelobe** artifacts
- Simplified imaging workflow for **a single SR mode** without the need to disjointly change between wide-field mode and SR mode as in HiFi-SIM
- Be suitable to combine with better parameter estimation algorithms to improve the accuracy of reconstructed images (e.g. PCA, Deep Learning)
- It's anticipated to even increase the overall speed by reducing the number of frames required to produce an SR image from nine to four

References

- [1] Rapid, artifact-reduced, image reconstruction for super-resolution structured illumination microscopy
<https://www.the-innovation.org/article/doi/10.1016/j.xinn.2023.100425>
- [2] A Sidelobe Suppressing Beamformer for Coherent Plane Wave Compounding
<https://www.mdpi.com/2076-3417/6/11/359>
- [3] Argo-SIM
<https://www.klv.co.jp/product/slide-microscope/argo-hm-v2.html>
- [4] Superresolution structured illumination microscopy reconstruction algorithms: a review
<https://www.nature.com/articles/s41377-023-01204-4>
- [5] Wiener Filter
<https://www.youtube.com/watch?v=Nm0UXL1OpRE>
- [6] High-fidelity structured illumination microscopy by point-spread-function engineering
<https://www.nature.com/articles/s41377-021-00513-w#Sec15>
- [7] Blind Deconvolution with Richardson Lucy Algorithm
<https://www.youtube.com/watch?v=v1kE3mbexko>
- [8] Scikit-image: image deconvolution
https://scikit-image.org/docs/0.23.x/auto_examples/filters/plot_deconvolution.html